



Vertical Sterilizer

(Pulsating vacuum)



Features

- ☆ Hand wheel translation and quick opening structure
- ☆ Three times pre-vacuum
- ☆ High quality stainless steel material
- ☆ Microcomputer automatic control, arbitrary setting sterilization parameters
- ☆ Safety interlocking device
- ☆ LED touch screen display
- ☆ Standard test interface
- ☆ Self-expanding seal
- ☆ With drying function
- ☆ With automatic protection function: over temperature protection; over pressure self discharge protection ; low water level protection; anti dry burning.
- ☆ With stainless steel bucket or basket
- ☆ Buzzer reminder after sterilization, automatic stop.
- ☆ Automatic discharge of cold air, automatic steam exhaust after sterilization.
- ☆ Internal printer



Model	YJ-MD35	YJ-MD50	YJ-MD75	YJ-MD100
Technical parameter				
Volume	35L	50L	75L	100L
Power	3.5kW	3.5kW	4.5kW	4.5kW
Voltage	AC 220V、50HZ			
Rated working pressure	0.22Mpa			
Rated working temperature	134℃			
Sterilization temperature	50℃~134℃			
Sterilization time setting range	4~120 min			
Drying time setting range	0~240min			
Chamber volume	Ø350×400	Ø400×505	Ø400×655	Ø450×650
Bucket dimension	Ø330×320	Ø380×430	Ø380×560	Ø420×540
Basket dimension (optional)	Ø320×350	Ø360×280×1	Ø360×280×2	Ø410×300×2
Outer dimension	698×498×940	630×530×1180	630×530×1180	760×560×1260
Packing size (L× W× H)mm	790×610×1130	800×600×1230	840×610×1310	870×650×1380
Gross/net weight	130/115KG	138/118KG	150/124KG	159/132KG